


In touch with Rheumatology



MyLabOne[™]
Musculoskeletal



- Dedicated ergonomics
- High-resolution imaging
- Simple and intuitive touch screen interface
- High Sensitivity Color and Power Doppler
- On-board MyLibrary reference guide
- Dedicated accessories and tools
- Up-to-date connectivity, including wireless

- Up to 22 MHz Transmission Frequency
- **XView** Image Enhancement
- **MView** Spatial Multi-View Imaging
- Hi-Sensitivity Power Doppler
- Measurements under 1 mm with Two Decimal Digits

A universe under the millimeter™




MyLibrary


Made in cooperation with Prof. Walter Grassi - Cattedra di Reumatologia, Università Politecnica delle Marche, Ancona, Italy.

Probes range


Remote controls integrated on the probe (customizable)



Flexor Pollicis Longus Tendon



Wrist, volar longitudinal



Joint inflammation, high-sensitivity Power Doppler

169 0043 00 (MA Rev.01)